

Ear Wash System

An easy-to-afford, easy-to-use ear irrigation system to remove cerumen without the mess.

Ear washing is one of the most commonly performed procedures in the primary care office. The Welch Allyn Ear Wash System provides a safe and effective way to remove cerumen from patients of all ages. This method of irrigation is much cleaner and safer than other irrigation methods currently on the market. Since it is easy to use with built-in safety and simplicity, it is able to be used by any healthcare professional within your facility.

Features

- Simplicity of operation — quick and easy set up.
- Pressure regulation — regulates the output pressure to prevent damage to the tympanic membrane or ossicular chain.
- Temperature indicator/control — light will turn green when the water is at body temperature, and will automatically turn off if the water gets too hot.
- Suctioning action to reduce messy “outflow”.
- Multi-functional disposable tip provides a seal at the ear canal. It has an outflow channel to direct water back to the drain, and a filter in the tip collects the cerumen, allowing the practitioner to see when the cleaning is complete. The tips can be used on children and adults.
- Portability — can be moved from room to room easily or transported off-site if necessary. All necessary components store easily in the storage tray.
- Affordable.

Model Number

29300*: Welch Allyn Ear Wash System

Accessory Numbers

- 29360: Eartips for Ear Wash System (case of 4 boxes or 100 eartips)
- 29330: Ear Wash System — Tubing
- 29371: Aerator Adapter Kit (includes snap aerator and 3 adapters)
- 29355: International Aerator Adapter Kit

Warranty

1 year warranty.

Patents

Patent pending.

*For appropriate product number outside of U.S.A., please advise language preference and country. Languages offered are English, French, Italian, German, Spanish, Portuguese, Chinese, and Dutch.



Best Ear Cleaning Solution



Easy To Use